

## Bath fluid Thermal P60

JULABO Thermal bath fluids are ideally suited for all of your temperature control applications and guarantee safe and reliable operation. Choosing the proper bath fluid is critical for best results in temperature control. Viscosity, oxidation characteristics and heat transfer of Thermal fluids are specifically matched with each JULABO temperature control unit.

JULABO Thermal bath fluids have been carefully selected and are optimized for maximum performance with JULABO instruments.

### Benefits

- Low toxicity
- Low viscosity
- High stability
- Minimum odor
- Good heat conductivity
- Low corrosion tendency
- Wide temperature ranges



### Technical Data

|                                         |            |
|-----------------------------------------|------------|
| Order no. 5 liters                      | 8891413    |
| Order no. 10 liters                     | 8891412    |
| Working temperature (°C)                | -60...+250 |
| Flash point (°C)                        | +130       |
| Fire point (°C)                         | +148       |
| Viscosity (kinematic at 20 °C in mm²/s) | 5          |
| Density (at 20 °C in g/cm³)             | 0.918      |
| Pour point (°C)                         | -65        |
| Boiling point (°C)                      | >205       |
| Ignition temperature (°C)               | 350        |
| Color                                   | clear      |