CORIO CD-BC6 Heating circulator



for internal and external temperature applications with stainless steel bath tanks and working temperatures from $+20~^{\circ}\text{C}$ to $+150~^{\circ}\text{C}$

Product features

- Models for internal and external applications
- Bright, white, easy to read display
- Very quiet
- Easy switching between internal and external circulation
- USB connection
- High-quality bath tanks made of stainless steel with bath lid and drain tap
- External pump connections
- Integrated cooling coil for counter-cooling
- Easy pump change-over between internal and external circulation
- For internal and external applications
- Precise temperature control



115V/60Hz (Nema N5-15 Plug)		
Heating capacity kW	1	
Viscosity max. cST	50	
Pump capacity flow rate I/min	16	
Pump capacity flow pressure psi	4.8	
Power A	10	

230V/60Hz (Schuko Plug - CEE 7/4 Plug Type F)	
Heating capacity kW	2
Viscosity max. cST	50
Pump capacity flow rate I/min	17
Pump capacity flow pressure psi	6.2
Power A	10

Technical data

Available voltage ver	sions	Bath	
Order No.	9 012 506	Bath tank	Stainless steel
Available voltage versions	x:	Bath cover	integrated
9012506.01	100V/50-60Hz (Nema N5-15 Plug)	Usable bath opening in. (W x L / D)	5.1 x 5.9 / 7.9
9012506.02	115V/60Hz (Nema N5-15 Plug)		
9012506.03	230V/50Hz (Schuko Plug - CEE 7/4 Plug Type F)		
9012506.03.chn	230V/50Hz (CN Plug)		





Product data sheet - CORIO CD-BC6



9012506.04	230V/50Hz (UK Plug Type BS1363A)
9012506.05	230V/50Hz (CH Plug Type SEV 1011)
9012506.13	230V/60Hz (Schuko Plug - CEE 7/4 Plug Type F)

Cooling		Other	Other	
Cooling coil integrated	Classification	Classification III (FL)		
	IP Code	IP 21		
	Pump function	Pressure Pump		
		Pump type	Immersion Pump	

Electronics		Dimensions and volumes	
Digital interface	USB	Weight lbs	22
Temperature control	PID1	Barbed fittings inner diameter	8/12 mm
Absolute temperature calibration	1 Point Calibration	Dimensions in. $(W \times L \times H)$	9.4 x 17.3 x 18.5
Temperature displayTemperature display	LED	Filling volume I	4.5 6
Temperature settingTemperature setting	Keypad	Pump connections	M16x1 male

Temperature values		Included in delivery
Setting the resolution of the temperature display °C	0.1	2 Barbed fittings for tubing 8 and 12 mm ID. (Pump connections M16x1 male)
Working temperature range °C	+20 +150	
Temperature stability °C	±0.03	
Ambient temperature °C	+5 +40	
Temperature display resolution °C	0.01 0.1	

For applications near or below ambient temperature: use a cooling coil or JULABO immersion cooler.

All Benefits



100% Checked.

100% testing. 100% quality. Each JULABO Circulator undergoes thorough quality testing before leaving the factory.



ATC.

Absolute Temperature Calibration, 1-point calibration (CD).



Handle with ease.

Makes day-to-day work easy. Comfortably move your CORIO around by using the ergonomic handles (front and rear).



Condensation protection.

Superb design solution. Integrated ventilation directs air over the bath lid and minimizes condensation.



Green technology.

Development consistently applied environmentally friendly materials and technologies.



Internal. External.

The pump is controlled via a lever located directly below the display. Easily change between internal and external circulation.



JULABO. Quality.

Highest standards of quality for a long product life.



Modern. Reliable.

High-grade components in every JULABO Circulator – platinum sensors, proven motor technology, CAN-Bus communication and much more.







Solid.

Minimized energy loss through high-quality insulation.



Stable.

Rubber feet allow for a secured footing of your CORIO to prevent damage to your laboratory equipment.



Tidy.

The special drain tap for easy draining of bath fluids without tools.



Touching permitted.

Optimum safety. The ergonomic plastic handle protects your fingers from hot surfaces.



Quick start.

Individual JULABO consultation and comprehensive manuals at your disposal.



Timer. Integrated.

CORIO circulators include an integrated timer function. When the set time has elapsed, a signal sounds and the device switches off. Setting range: 0 ... 999 minutes.



Brilliant.

Very bright display makes it easy to read even from a distance.



Connection. Easy.

Inclined pump connections (M16×1) facilitate the connection of applications. Each unit includes 2 barbed fittings of 8/12 mm diameter each.



Everything at the front.

All operating controls and safety functions are accessed easily and comfortably from the front.



Exact.

You can rely on it. PID1 control and 'Active Cooling Control' make the new CORIO precise and perfect.



Locked in.

The lockable power plug guarantees safe connection. More process safety.



Switch on. And off you go.

Intelligent operating concept. Ready for operation with just a few quick and easy steps.



Early warning system for low liquid level

Maximum safety for applications, optical and audible alarm, allows user to refill bath fluid before the unit shuts down



Connectivity.

Remote control made easy. CORIO CD circulators feature a USB connection.

