QUICK. CONFIDENT. RELIABLE. DOWN TO THE CORE.

It's simple. Quickly remove the valve core and reduce the time it takes to recover refrigerant and pull a system vacuum. The hilmor Valve Core Removal Tool makes quick work of removing and replacing Schrader valves. A unique, spring gripper securely grabs the valve core every time and increases the durability of the tool. And the built-in plunger magnet allows for convenient storage once the core is removed from the system. Because when you're pulling a vacuum and recovering refrigerant, nothing should get in the way of your flow.

Valve Core Removal Tool Instructions:

- Pull back plunger and connect the female fitting of the Valve Core Removal Tool to the system access port.
- 2. Rotate the plunger while pushing toward the valve core until the gripper securely engages. Push firmly to grip the valve core.
- 3. Rotate the plunger counter-clockwise to unscrew the valve core from the system.
- 4. Pull plunger back until it is fully extended. (It is common for the system pressure to push the plunger toward you.)
- 5. Close the ball valve handle.
- Unscrew the plunger from the body of the tool to remove the plunger with the valve core.
- 7. Connect your refrigeration hose to the Valve Core Removal Tool body.
- 8. Open the ball valve handle to get the job done.

Note:

When evacuating a system, connect a vacuum gauge to the side port.

hilmer

