









# Precision Monitoring with NIST Traceability and Compliance Built-In

EDGE Pharma is engineered for strict compliance in pharmaceutical, vaccine, and cold chain applications. Supplied with NIST-traceable external temperature probes, it ensures reliable, audit-ready accuracy and meets NSF/ANSI 456 standards.

The device clearly displays NIST certificate validity onscreen and sends automatic email reminders before expiration, keeping your monitoring system compliant with minimal manual oversight.

With seamless connectivity, powerful cloud features, and smart alarm handling, EDGE Pharma offers a complete, worry-free solution for high-stakes temperature monitoring.

NIST validity: 3 years in use (5 years on shelf), with automatic reminders shown on-screen and sent via email. Renew certification easily by replacing the affordable probe, no need to replace the entire device.







Cellular

**Sensors & Monitoring** 



2 x Temp Probes + NIST



**Ambient** Temperature



Relative Humidity



Power Disruption



Digital Input



Mold Index



Location Cellular Only

- Integrated Ambient **Temperature** and **Relative Humidity** Sensor.
- Two (2) External Digital Temperature Probe Inputs (1.5m, 5ft.).
- NIST-Traceable Calibration Certificate for External Probes.
- One (1) **Digital Input** (Dry Contact) for Door Magnet or Flood Sensor.
- Up to 2 years **battery life** (4 x AA replaceable Alkaline)

#### **Real-time Alarms**

- **High & Low Temperature** and Humidity alarms for All Four Channels.
- **Power Interruption** and Restoration Alarms.
- Flood Detection, Door Magnet or Dry contact Real-Time Alarms.
- Mold Index Calculation and Mold risk Warnings.









Audible

Visible

LOGGERFLEX

Notice of Manifestories Solutions





### **EDGE's Special Features**

#### Advanced Monitoring, Alarms, and Seamless Connectivity

Connects directly to WiFi or Cellular (No gateway required) Alarms via Email, SMS, Phone Call

Power Disruption and Re-connection Notification NSF/ANSI 456 compliant Audible and Visual alarms

Supports NIST-traceable calibration with automatic expiration alerts

Real-time BLE streaming

Automatic WiFi-to-cellular rollover (in Cellular model)

Works on 4xAA batteries (up to 2 years) or USB-C power

Super user-friendly setup with QR scan

#### Display & Interface

Graphical e-paper screen shows:

- Real-time data
- Active alarms
- Battery, Connection and power status
- NIST certificate expiration

Audible and visual alarms directly on device.

## **LoggerFlex Cloud Platform**

## Powerful Web App with Compliance-Built Features

Automatic Temperature Logbook
FDA CFR21 Part 11 Compliance
Secure, untamperable activity log
Time-stamped user activity and alarm acknowledgments
Professional Reporting with:

- Min/Max records
- MKT (Mean Kinetic Temperature)
- Graphs and analytics
- Geographical Management (Map View)
   Dedicated Mobile & Desktop Interfaces
   Team Sharing & Role-Based Access Control

#### Connectivity

WiFi for continuous cloud sync and Real-tima alarms BLE for smooth real-time data broadcast Cellular with automatic WiFi failover. (Cellular model)

Designed to meet the performance requirements of the **NSF / ANSI 456** Vaccine Storage Standard.

# **Technical Specifications**

I and the second		
Built in sensors		Temperature and Relative Humidity (RH)
Operation Temperature range	°C	-20 to +70
	°F	-4 to +160
Temperature Measurement Accuracy		±0.2 °C (±0.36 °F)
Temperature Reporting Resolution		0.1
RH measurement range		0-99% non-condensing
RH Measurement Accuracy		±1.8% RH (at 25 °C, 10–90%)
Interface		Wi-Fi - IEEE 802.11 b/g/n – 2.4 GHz
		BLE (Bluetooth Low Energy) – 2.4 GHz
		Cellular LTE-M & 2G (all band - Global) *
Temperature Probes (x2) Measurement Range	°C	-55 to +125
	°F	-67 to +257
Temperature Probes (x2) Measurement Accuracy		±0.5 °C (±0.9 °F)
FCC ID	WiFi	2AC7Z-ESPWROOM32
	Cellular	2AJYU-8VC0001
Power Supply	Internal	4 x AA batteries
	External	5V DC Standard USB-Charger
Display Update Temperature Range Out of this range, display updates are paused	°C	0 to +50
	°F	+32 to +122
Internal Memory Capacity		49,000 Record of each measured Parameter
Recoding interval		1 to 30 minutes
Syncing interval (sending data to the cloud)		Real-time Alarms - 10 minutess to 7 days
		A

<sup>\*</sup> Cellular model - with global e-sim





