

# P3000 Series

## Pneumatic Deadweight Testers Model P3000

### Technical Data

#### Features

- Pressure ranges from Vacuum through 2 000 psi (140 bar)
- Accuracy better than 0.015 % of reading. (Increased accuracy option of 0.008 %)
- Gas operated with optional hand pumps
- Dual vacuum/pressure models available
- Psi, bar, kgf/cm<sup>2</sup>, kPa and MPa ranges available
- Piston/cylinder design provides stability and repeatability
- Built-in pneumatic hand pumps for pressure and vacuum
- High quality needle valves provide optimum control
- Mounted spirit level and adjustable feet
- Test station design with O-rings eliminates the need for PTFE tape or wrenches
- Improved case with spring-loaded latches for hood
- Sturdy weight box with hinged lid and side handles for easy transport



The P3000 Series is the culmination of over 50 years experience in the production and design of primary pressure standards. With features designed to improve accuracy and performance, increase reliability and simplify operation these deadweight testers can be used to calibrate virtually any pressure sensing device, including transducers, transmitters, gauges or pressure switches.

Each instrument is supplied with a detachable lid that makes it neat, compact and easily portable. The masses are stored in a high-quality case with

a self-locking mechanism to protect them during transit.

All units are provided with an accredited calibration certificate, M14X1.5, M20X1.5, 1/8, 1/4, 3/8, and 1/2 NPT and BSP female adaptors, operating fluid (where applicable) and spare seals.

The piston/cylinder assemblies are manufactured to the very highest standards and are supplied with accredited calibration reports with stated instrumental measurement uncertainties.

## Operating principle

Deadweight testers are the primary standard for pressure measurements. Utilizing the well proven piston-gauge system consisting of a vertically mounted, precision lapped piston and cylinder assembly, accurately calibrated masses (Force) are loaded on the piston (Area), which rises freely within its cylinder. These masses balance the upward force created by the pressure within the system.

## Instrument base


There are four basic variations in the pneumatic offering; single pressure, single vacuum, dual vacuum and pressure, and oil-lubricated pneumatically driven high-pressure models. Pressure or vacuum is generated by an external source or an optional built-in hand pump (standard on the lower pressure models). The hand pump option is not available on instruments above 500 psi (3.5 MPa).

## Piston/cylinder assemblies

The piston/cylinder assembly is the heart of each deadweight tester. They are manufactured from materials that provide stability, durability, and low thermal coefficients and distortion. Our experience and knowledge of piston/cylinder production and calibration ensure the precision and performance required for today's demanding calibration requirements.

## Masses

Standard masses are series 3 non-magnetic austenitic stainless steel. Each mass is marked with the serial number of the instrument and the nominal pressure value relative to the high- or low-pressure piston, when applicable. The masses for the vacuum ranges and the optional fine increment masses are stainless steel and/or aluminum.



## Gravity correction

Gravity varies significantly with geographical location and this variation has a direct effect on the force of the weights and the accuracy of the deadweight tester. Each instrument can be calibrated to local gravity at no extra cost. If unspecified, instruments will be calibrated to Standard Gravity at 980.665 cm/s<sup>2</sup>.

## Gas supply

All instruments are fitted with a 1/4 inch NPT female external supply port for connection to the gas supply source. A compressed gas bottle (nitrogen or dry air) fitted with a pressure regulator is recommended. Ten percent above the maximum pressure is required. A built-in hand pump is available on some models to generate pressures to 300 psi (2 MPa).

For vacuum requirements, an external vacuum pump can be connected to the 1/4 inch port. The built-in hand pump will also generate up to 90% vacuum.


## PressCal Software

Windows-based software program that allows users to easily apply all necessary corrections to enhance the deadweight tester performance. Calibration details are then stored and/or used to automatically create a calibration certificate.

The use of PressCal allows for an improved accuracy of 0.008%.

### P3000 Series Vacuum and Dual Vacuum/Pressure Models

The combined vacuum/pressure models are extremely versatile instruments providing vacuum through 500 psi (3.5 MPa) capacity in one unit. An optional built-in hand pump can be installed to generate vacuum (90 %) through 300 psi (2 MPa).




### P3000 Series Liquid Lubricated Gas Operated Models

The Models P3031 and P3032 feature an oil lubricated piston/cylinder assembly, which ensures that the instrument is less susceptible to performance issues caused by contamination of the gas supply or the environment in which the unit is being operated.

The supply gas, via the deadweight tester controls, is introduced into the piston/cylinder chamber. As the gas acts directly on the surface of the oil, there is a rate-less transfer of pressure within the system.

The instrument is designed to prevent inadvertent over-filling of the chamber, thus avoiding contamination of the gas system. There is no physical barrier between the gas and oil so there is some "oil-mist" transfer potential. Therefore, we do not recommend the standard version for oxygen service. A special version of these instruments using a hydro-carbon free fluid is available.



Cross-sectional illustration of liquid lubricated piston chamber



## Specifications

| Pressure ranges                                                            |                                                                                                                                                                                                                                                        |                   |                          |
|----------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------|--------------------------|
| Instrumental measurement uncertainty                                       | ± 0.015 % of reading (± 0.008 % optional)<br>Note: For more information on the uncertainty, including the % reading threshold, please see technical note, "Guide for the Uncertainty Analysis in Pressure When Using P3000 Series Dead-weight Testers" |                   |                          |
| Mass and PCU materials of construction                                     |                                                                                                                                                                                                                                                        |                   |                          |
| Standard weight material                                                   | Series 3 non-magnetic austenitic stainless steel<br>Density: 7.8 g/cm³                                                                                                                                                                                 |                   |                          |
| Optional fine increment masses                                             | Aluminum<br>Density: 2.7 g/cm³                                                                                                                                                                                                                         |                   |                          |
|                                                                            | Piston material                                                                                                                                                                                                                                        | Cylinder material | Coefficient of expansion |
| P3011, P3012, P3013, P3022, P3023, P3025 (V)                               | Ceramic                                                                                                                                                                                                                                                | Martensitic steel | 11 ppm/°C                |
| P3014, P3015, P3025 (P), P3031, P3032 (P & V)                              | Tungsten carbide                                                                                                                                                                                                                                       | Martensitic steel | 16.5 ppm/°C              |
| General                                                                    |                                                                                                                                                                                                                                                        |                   |                          |
| Test port adaptors                                                         | M14 X 1.5, M20 X 1.5, 1/8 in, 1/4 in, 3/8 in and 1/2 in NPT and BSP                                                                                                                                                                                    |                   |                          |
| Weight                                                                     | 11 kg (24 lb)                                                                                                                                                                                                                                          |                   |                          |
| Dimensions (W x D x H)                                                     | 440 mm x 300 mm x 215 mm (17.5 in x 12 in x 8.5 in)                                                                                                                                                                                                    |                   |                          |
| Mass set weight (typical)                                                  | 29 kg (65 lb)                                                                                                                                                                                                                                          |                   |                          |
| Weight increments                                                          |                                                                                                                                                                                                                                                        |                   |                          |
| Minimum standard weight increments                                         |                                                                                                                                                                                                                                                        |                   |                          |
| P3012, P3022 Pressure                                                      | 1 inH <sub>2</sub> O, (5 mbar, 0.5 kPa)                                                                                                                                                                                                                |                   |                          |
| P3013, P3023 Pressure                                                      | 2 inH <sub>2</sub> O, (10 mbar, 1 kPa)                                                                                                                                                                                                                 |                   |                          |
| P3014, P3015, P3025 Pressure                                               | 1 psi, (0.1bar, 10 kPa)                                                                                                                                                                                                                                |                   |                          |
| P3031, P3032                                                               | 1 psi, (0.1bar, 10 kPa)                                                                                                                                                                                                                                |                   |                          |
| P3011, P3022, P3023, P3025 Vacuum                                          | 0.2 inHg (10 mbar, 1 kPa)                                                                                                                                                                                                                              |                   |                          |
| Optional fractional weights                                                |                                                                                                                                                                                                                                                        |                   |                          |
| P3014, P3015, P3025 Pressure                                               | 0.1 psi (10 mbar, 1 kPa)                                                                                                                                                                                                                               |                   |                          |
| Built-in hand pumps                                                        |                                                                                                                                                                                                                                                        |                   |                          |
| Pressure mode                                                              | Max pressure 300 psi (2 MPa)                                                                                                                                                                                                                           |                   |                          |
| Vacuum mode                                                                | To 90 % vacuum                                                                                                                                                                                                                                         |                   |                          |
| Operating fluid for liquid lubricated models                               |                                                                                                                                                                                                                                                        |                   |                          |
| 55-65S mineral oil. For oxygen clean applications, please consult factory. |                                                                                                                                                                                                                                                        |                   |                          |

## Ordering information

### Air operated - single PCU

#### Model

**P3012-KPA-P** 1.5 to 100 kPa  
**P3012-MBAR-P** 15 to 1 000 mbar  
**P3012-INH2O-P** 5 to 400 inH<sub>2</sub>O  
**P3013-KPA-P** 3 to 200 kPa  
**P3013-MBAR-P** 30 to 2 000 mbar  
**P3013-INH2O-P** 12 to 800 inH<sub>2</sub>O  
**P3014-BAR-P** 0.2 to 10 bar  
**P3014-KGCM2-P** 0.2 to 10 kgf/cm<sup>2</sup>  
**P3014-PSI-P** 3 to 150 psi  
**P3014-KPA-P** 20 to 1 000 kPa  
**P3014-5MPA-P** 0.02 to 1 MPa  
**P3015-BAR-P\*** 0.2 to 35 bar  
**P3015-KGCM2-P\*** 0.2 to 35 kgf/cm<sup>2</sup>  
**P3015-PSI-P\*** 3 to 500 psi  
**P3015-KPA-P\*** 20 to 3 500 kPa  
**P3015-MPA-P\*** 0.02 to 3.5 MPa

\* -P indicates built-in pump optional.

### Air operated - dual PCU

| Model                 | Pressure range                | Vacuum range     |
|-----------------------|-------------------------------|------------------|
| <b>P3022-KPA-P*</b>   | 1.5 to 100 kPa                | 3 to 100 kPa     |
| <b>P3022-MBAR-P*</b>  | 15 to 1,000 mbar              | 30 to 1 000 mbar |
| <b>P3022-INH2O-P*</b> | 5 to 400 inH <sub>2</sub> O   | 1 to 30 inHg     |
| <b>P3023-KPA-P*</b>   | 3 to 200 kPa                  | 3 to 100 kPa     |
| <b>P3023-MBAR-P*</b>  | 30 to 2 000 mbar              | 30 to 1 000 mbar |
| <b>P3023-INH2O-P*</b> | 12 to 800 inH <sub>2</sub> O  | 1 to 30 inHg     |
| <b>P3025-BAR-P*</b>   | 0.2 to 35 bar                 | 30 to 1 000 mbar |
| <b>P3025-KGCM2-P*</b> | 0.2 to 35 kgf/cm <sup>2</sup> | 30 to 760 mmHg   |
| <b>P3025-PSI-P*</b>   | 3 to 500 psi                  | 1 to 30 inHg     |
| <b>P3025-KPA-P*</b>   | 0.02 to 3.5 MPa               | 3 to 100 kPa     |
| <b>P3025-MPA-P*</b>   | 0.02 to 3.5 MPa               | 3 to 100 kPa     |

\* -P indicates built-in pump optional.

### Options

|                                                                            |
|----------------------------------------------------------------------------|
| Built-in pump (where applicable)                                           |
| Fine Increment Weights (where applicable)                                  |
| Conversion Weight Set to allow pressure increments in other pressure units |
| 0.008 % Improved Accuracy (PressCal Software)                              |
| COMPASS for Pressure Calibration Software                                  |

### Liquid lubricated

#### Model

**P3031-BAR** 1 to 70 bar  
**P3031-KGCM2** 1 to 70 kgf/cm<sup>2</sup>  
**P3031-PSI** 10 to 1 000 psi  
**P3031-KPA** 100 to 7 000 kPa  
**P3031-MPA** 0.1 to 7 MPa  
**P3032-BAR** 1 to 140 bar  
**P3032-2KGCM2** 1 to 140 kgf/cm<sup>2</sup>  
**P3032-PSI** 10 to 2 000 psi  
**P3032-KPA** 100 to 14 000 kPa  
**P3032-MPA** 0.1 to 14 MPa

### Vacuum operated - single PCU

#### Model

**P3011-KPA-P\*** 3 to 100 kPa  
**P3011-MBAR-P\*** 30 to 1 000 mbar  
**P3011-INHG-P\*** 1 to 30 inHg  
**P3011-MMHG-P\*** 30 to 760 mmHg

\* -P indicates built-in pump optional.

## Fluke Calibration.

*Precision, performance, confidence.™*

|                 |
|-----------------|
| Electrical      |
| RF              |
| Temperature     |
| <b>Pressure</b> |
| Flow            |
| Software        |

©2011 Fluke Calibration.  
 Specifications subject to change without notice.  
 Printed in U.S.A. 7/2012 3833551D D-EN-N

Modification of this document is not permitted  
 without written permission from Fluke Calibration.